

# NOCOLYSE+

The ready-to-use bio-disinfectant.



Manufactured in France by Oxy'Pharm according to ISO 13485.  
Use biocidal products carefully. Read the label and product information before use.

**Nocolyse+** is a surface bio-disinfection product. It is a ready-to-use solution with a base of 12% hydrogen peroxide, which must be used with spray devices from the **Nocotech** range. The use of both **Nocolyse+/Nocospray** (or **Nocomax**) is effective against all types of micro-organisms: this pairing disinfects surfaces with a bactericidal, fungicidal, virucidal, yeasticidal, tuberculocidal and sporicidal effect.



Biodegradable



No residue



Non-toxic



Non-corrosive



Non-allergenic



No germ resistance

## COMPOSITION

Stabilised hydrogen peroxide solution 12% (120 ml/l)  
• EC=231-765-0 / CAS=7722-84-1. Silver 17 ppm •  
EC=231-131-3 / CAS=7440-22-4.

## STORAGE

Store the product in its original packaging, upright and in a cool, well-ventilated area.

- Storage in the original closed packaging: 2 years from the date of manufacture.
- Storage once opened: 2 months from the date of opening.

## PRECAUTIONS FOR USE

Consult the product's safety data sheet, available on request by email: [info@oxypharm.net](mailto:info@oxypharm.net)

## REFERENCES AND PACKAGING

	NEUTRAL	MINT	NOCODOR
1 L BOTTLE	4040.001	4041.001	4042.001
6 X 1 L BOX	4040.001-6	4041.001-6	4042.001-6
5L CAN	4040.005	4041.005	4042.005
10L CAN	4040.010	4041.010	4042.010
20L CAN	4040.020	4041.020	4042.020



## INSTRUCTIONS FOR USE

- 1 Follow the instructions for use of the **Nocotech** spray device (see instructions for use and quickstart document).
- 2 Screw the 1 L bottle into the **Nocospray** or the 20 L can into the **Nocomax**.
- 3 Set the volume (V) according to the desired treatment.
- 4 After the end of diffusion, respect the minimum contact time indicated in the efficiency table below.

### IMPORTANT

- Throughout the operation of the machine and the contact time, keep the room closed and do not enter it. The treatment must be performed away from humans.
- Obtaining a good quality of disinfection is directly linked to the respect of a strict cleaning protocol, performed before the treatment.