



Manufactured in France by Sanivap according to ISO 13485. SANIVAP devices are Class IIa medical devices - Notified body: GMED 0459

Θ



SV10 - LUBRICANT

The **SV10** is a tube of FDA silicone lubricant. It gives proper connection for accessories.

COMPATIBILITY

• Steam hose ring O seals.

FIND OUT MORE

• Apply the product to clean seals.

05/16

-T-Consumable

Avoid excessive lubrication.

PRODUCT
Lubricant

REFERENCES FRANCE/WORLD SV10/8650.10



SV20 - LIMESTONE SEQUESTERING AGENT

SV20 is a limescale sequestering agent for boiler feed water (1L canister) to be diluted according to the dosage indicated on the label.

PRODUCT	REFERENCES FRANCE/WORLD
Limestone Sequestering Agent	SV20/8650.20

SP34 - SEAL KIT



The **SP34** seal kit has been designed for the maintanance of the appliance's seals (once a year).

CONSUMABLES Seal kit

REFERENCES FRANCE/WORLD SP34/8650.34

SV75 - FILTRATION CARTRIDGE

The SV75 device is a replacement cartridge for the SV310 filtration system.

COMPATIBILITY

• SV310 filtration system.

FIND OUT MORE

• It needs replacing once a year or when the volume of water consumed is reached.

PRODUCT Filtration cartridge REFERENCES FRANCE/WORLD SV75/8640.75

MV100 - ACTIVE CLEANER

MV100 is an active surface cleaner. It removes biofilm and encrusted residues from hard surfaces.

FIND OUT MORE

- Comes with a spray.
- Do not use on aluminium, fabric or fragile surfaces.
- Packaging: 500mL canister.



PRODUCT

Active Cleaner

REFERENCES FRANCE/WORLD

MV100/8650.100





MV200 - ACTIVE FOOD CLEANER

MV200 is an active surface cleaner that can come into contact with food. It removes biofilm and encrusted residues from hard surfaces.

FIND OUT MORE

- Comes with a spray.
- Do not use on aluminium, fabric or fragile surfaces.
- Packaging: 500mL canister.

PRODUCT

Active food cleaner

REFERENCES FRANCE/WORLD

MV200/8650.200



MV600H-2L

MV600H is a detergent solution that can be used for applications in contact with food . The MV600H product must only be used with the SP540H device.

Composition: Stabilised hydrogen peroxide 7.9% - EC = 231-765-0 / CAS = 7722-84-1

The solution is stable for 24 months in the original closed packaging from the manufacturing date. Stable from 2 months from the date of opening the packaging.

PRODUCT

H₂O₂-based solution

REFERENCES FRANCE/WORLD

MV600H-2L/8670.600H-2L

